Name:	
Signature	Centre /Index No

545/3 Chemistry Paper 3 July/August,2022 2 Hours.



TESO SECONDARY SCHOOLS MOCK EXAMINATION ASSOCIATION(TESSMEA)

Uganda Certificate of Education 2022

CHEMISTRY PRACTICAL)

Paper 3

2 Hours.

Instructions to candidates

Answer both questions. Answers are to be written in the spaces provided in this booklet

All your working must be in blue or black ink. Any work done in pencil will not be marked

You are not allowed to use any reference books (i.e text books, booklets on qualitative analysis etc).

All working must be clearly shown.

For Ex	aminer's use only	18 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1
Q.1		
Q.2		
Total		

Q1.	You are provided with the following BA1, which is a solution made by dissolving 1g of sodium hydroxide in 250cm ³ of water.			
	BA2, which is a solution containing 6.3g dm ⁻³ of an acid H ₂ A .x H ₂ O.			
	You are required to determine the value of X in the acid			
	Procedure:			
	Pipette 25cm ³ (or 20cm ³) of BA2 into a conical flask. Add 2-3 drops of phenolphthalein indicator and titrate with BA1 from the burette.			
	Repeat the titration until you obtain consistent results			
	Record your results in the table below			
	Volume of pipette used			
	Final burette reading (cm ³)			
	Initial burette reading(cm ³)			
	Volume of BA1 used (cm ³)			
	Titre values of BA1 used for average volume. (01 mark) cm ³)		
	Average volume of BA1 used. (2 ½ marks)		
	Questions			
	(a). Calculate			
	(i). The molarity of BA1 [Na = 23, H = 1, O = 16]. (03 marks)			

	(11).	The number of moles of BA1 that reacted with BA2. (1 $\frac{1}{2}$ marks)
		·
	(iii).	The number of moles of BA2 that reacted with BA1. (03marks)
	(iv).	The number of moles of the acid H ₂ A.x H ₂ O in 1.0dm ³ of BA2
		(1 ½ marks)
	(b).	Determine the relative formula mass of the acid $H_2A.xH_2O$, hence the value of X. [$H = 1$, $A = 88$, $O = 16$]. (4 ½ marks)
Q2.	anions	are provided with substance T, which contains two cations and two s. You are required to carry out the following tests on T and identify tions and anions in T. Identify any gas(es) evolved and record your vations and deductions in the spaces provided. (25marks)

Tests	Observations	Deductions
a). Heat a spatula endful of		
T strongly in a dry test tube		
	,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,	
h) To two enetule endfule		
b). To two spatula endfuls of T in a test tube add dilute		
nitric acid drop-wise until	• .	
there is no further change,	,	4
followed by dilute sodium		
hydroxide solution drop-		
wise until in excess. Filter		
and keep both the filtrate		
and the residue.		
c). To the filtrate add dilute	• , ,	
nitric acid drop-wise until	T. A.	
the solution is just acidic. Divide the solution into four		
parts.	•	
(i). To the first part of the		
acidic solution add dilute		
sodium hydroxide solution		
drop wise until in excess.	• •	
(ii). To the second part of		
the acidic solution add		at adjust the same
dilute ammonia solution drop-wise until in excess.		
drop-wise until in excess.		-1
(iii). To the third part of the	,	
acidic solution add sodium		-
sulphate solution.		
- ·		
(iv). Use the fourth part of	,	
the acidic solution to		
carryout a test of your own		- 4
choice to confirm the cation in the filtrate.		
in the intrate.		
	· ·	

d). Wash the residue and dissolve it in dilute nitric acid. Divide the solution into two parts		
(i). To the first part add dilute sodium hydroxide solution drop-wise until in excess		
(ii). To the second part of the solution add dilute ammonia solution drop- wise until in excess		
	andand	

END

545/3

CHEMISTRY

PRACTICAL

3 Hours



TESO SECONDARY SCHOOLS MOCK EXAMINATION ASSOCIATION (TESSMEA)

Uganda Certificate of Education 2022

CHEMISTRY PRACTICAL

Paper 3

CONFIDENTIAL

Great care should be taken that the information given below does not reach the candidates either directly or indirectly.

NB: The Headteacher must ensure that the teacher responsible for preparing the apparatus hands in his/her trial results properly sealed in a separate envelope and firmly fastened to the candidate's scripts envelope.

- 1 The description of the reagents and chemicals specified below does not necessarily correspond with the description in the question paper. Candidates must NOT be informed of the difference.
- 2 Candidates are not allowed to use reference books (i.e text books) booklets of qualitative analysis etc) during examination.
- 3 In addition to fittings and substances ordinarily contained in a chemistry laboratory, each candidate will require:
 - 1 burette (50cm³)
 - 1 pipette (25.0cm³ or 20.0cm³)
 - 1 conical flask
 - 1 measuring cylinder (50cm³)
 - 6 test tubes

1 filter paper

1 filter funnel

2 beakers (250ml)

100cm³ of distilled water

100cm³ of BA1

100cm³ of BA2

2 g of H

Easy access to:

- Heat source
- · Common reagents for identifying gases, cations, and anions
- Phenolphthalein and methyl orange indication

BA1 is prepared by dissolving 4.0g of sodium hydroxide in distilled water to make one litre of solution.

BA2 is prepared by diluting 8.6cm³ of concentrated hydrochloric acid ((36%, 1.18g cm³) with distilled water to make one litre of solution

H - mixture of magnesium carbonate and lead(II) nitrate in a ratio of 2:1